

CANEXIS™

Integrated Workflow Solution for Endoscopy

- Audit Preparedness
- Patient-Ready Endoscopes
- Efficient Workflow and Tracking



Are you prepared for a Joint Commission audit?



Are you confident that a dirty endoscope will never be used on a patient?



Is it a scramble to locate endoscopes and other equipment for upcoming procedures?

CANTEL HELPS YOU ANSWER THESE QUESTIONS AND MORE WITH THE CANEXIS™ INTEGRATED WORKFLOW SOLUTION.

The CANEXIS™ Integrated Workflow Solution was developed in close collaboration with endoscopy experts to create a user-friendly tool that gives you peace of mind that your facility is always executing best practices.

The software is flexible and configurable to support your unique workflow needs. Unparalleled quality management tools and the ability to achieve full traceability of your procedures and endoscopes make CANEXIS Integrated Workflow Solution the most useful tool on the market to manage your endoscopy department.





AUDIT PREPAREDNESS

- ✓ Are you prepared for your next Joint Commission audit?
- ✓ Can you demonstrate that your staff has completed all necessary steps in the reprocessing cycle?

Our system can help. It is designed to make sure your workflow follows all relevant societal guidelines and provides a record that every required step was followed during reprocessing.

The best part? It's fully automated. The CANEXIS™ Software communicates with AERs, drying cabinets and other equipment in your department for a complete picture of reprocessing. That means you know which patients an endoscope was used on, who reprocessed it, which AER it was reprocessed in and when it was last out for repair.

The reporting tool within the software allows you to immediately access the information you need to do your job or respond to an auditor's questions. The information is all there, right at your fingertips.

CANEXIS™ Integrated Workflow Solution for Endoscopy			
 Bedside	 Manual Clean	 Disinfection	 Storage
Notes the staff member who registered the patient, endoscope, physician and consumables	Confirms manual cleaning starts within one hour	Confirms the staff member who started or canceled the HLD cycle	Notes the staff member who loaded the endoscope for storage
Notes which staff member completed the bedside cleaning and begins timing to track Golden Hour compliance	Tracks leak testing, flushing and visual inspection and notes the outcome	Tracks whether the HLD cycle passed or failed	Tracks your hospital policy for hang/storage time

PATIENT-READY ENDOSCOPES

SOFTWARE TO SUPPORT YOUR WORKFLOW NEEDS

- ✓ Are you confident that a dirty endoscope will never be used on a patient?


As it relates to patient safety, this may be the most important question you have to consider.

Ensuring patient safety means ensuring patient-ready endoscopes are available for every procedure and for every patient.

But it can be difficult to know with certainty that all reprocessing steps are being completed within the correct time and in the correct order.

In the event of a missed reprocessing step, the software notifies your staff and prevents the endoscope from moving on in the workflow until the required step is performed. This information is recorded and can be used to develop training plans for your staff based on their performance.

ENDOSCOPE TRACKING EACH STEP OF THE WAY



EFFICIENT WORKFLOW AND TRACKING

✓ Can you locate your endoscopes?

✓ Can you meet your procedure needs each day?

The CANEXIS™ Integrated Workflow Endoscopy Solution helps you better manage your daily schedule and makes your whole department more efficient. At a glance, you can:

- Know which endoscopes are ready for use.
- Identify which endoscopes need to be reprocessed and where in your facility they are located.
- Track which endoscopes are out for repair.
- Monitor trends over time to help you make repair vs. replace decisions.
- Determine the status and location of your entire endoscope inventory.





THE COMPLETE CIRCLE OF PROTECTION

For more than 40 years Cantel has been delivering innovation and excellence in endoscope reprocessing. We offer a comprehensive line of endoscope reprocessing products including market-leading detergents and disinfectants, technologically advanced automated endoscope reprocessors, transportation systems, drying and storage cabinets and an integrated workflow solution.

Our ecosystem of product and service offerings comes together in our **Complete Circle of Protection**, an infection prevention program designed to help you streamline reprocessing workflow, improve department efficiency and reduce the risk of infection.

CANEXIS™ is a trademark of Cantel Medical Corp., its affiliates or related companies.

MEDIVATORS IS NOW CANTEL

www.medivators.com

TO PLACE AN ORDER

p: 1.800.328.3340 (customer service) f: 1.800.686.8493 e: custserv@medivators.com

© 2020 Cantel Medical Corp. LIT-00229/A

